

Item no	Item		
SatinSheet 2Direction			
IM4133S	SatinSheet 2Direction mini - small	W1950xL1050 mm (satin W600xL1050 mm). Stripes	€ 84,70
IM4108S	SatinSheet 2Direction mini	W1950xL1050 mm (satin W700xL1050 mm). Stripes	€ 73,57
IM4175S	SatinSheet 2Direction mini - wide	W2300xL1050 mm (satin W850xL1050 mm). Stripes	€ 99,99
IM4178S	SatinSheet 2Direction mini - double	W2300xL1050 mm (satin W960xL1050 mm). Stripes	€ 127,49
IM4119S	SatinSheet 2Direction midi - small	W2000xL1400 mm (satin W600xL1400 mm). Stripes	€ 97,67
IM4113S	SatinSheet 2Direction midi	W2000xL1400 mm (satin W700xL1400 mm). Stripes	€ 84,26
IM4176S	SatinSheet 2Direction midi - wide	W2300xL1400 mm (satin W850xL1400 mm). Stripes	€ 106,33
IM4179S	SatinSheet 2Direction midi - double	W2300xL1400 mm (satin W960xL1400 mm). Stripes	€ 136,05
IM4181S	SatinSheet 2Direction maxi - small	W2000xL2000 mm (satin W600xL2000 mm). Stripes	€ 156,57
IM4114S	SatinSheet 2Direction maxi	W2000xL2000 mm (satin W700xL2000 mm). Stripes	€ 106,82
IM4177S	SatinSheet 2Direction maxi - wide	W2300xL2000 mm (satin W850xL2000 mm). Stripes	€ 122,11
IM4180S	SatinSheet 2Direction maxi - double	W2300xL2000 mm (satin W960xL2000 mm). Stripes	€ 165,82
IM4122S	SatinSheet 2Direction 180 cm bed	W2800xL2000 mm (satin W600xL2000 mm). Satin not ce	€ 143,55
IM4112S	SatinSheet 2Direction Fit	W900xL2000 mm (satin W700xL2000 mm). Stripes	€ 119,45
IM4142S	SatinSheet 2Direction Elastic Fit	W900xL2000 mm (satin W700xL2000 mm). Stripes	€ 142,30
IM4167S	SatinSheet 2Direction Corner - small - mid	W900xL1730 mm (satin W670xL1730 mm). W/elastic cord	€ 129,23
IM4157S	SatinSheet 2Direction Corner - midi	W850xL2000 mm (satin W650xL1400 mm). W/elastic cord	€ 97,33
IM4107S	SatinSheet 2Direction Corner	W850xL2000 mm (satin W650xL2000 mm). W/elastic cord	€ 85,28
IM4147S	SatinSheet 2Direction Corner - wide	W1050xL2000mm (satin W800xL2000 mm). W/elastic cor	€ 136,58
IM4137S	SatinSheet 2Direction Corner - double	W1200xL2000mm (satin W900xL2000 mm). W/elastic cor	€ 118,48
IM4197S	SatinSheet 2Direction Corner - double - m	W1800xL2000mm (satin W1550xL2000 mm). W/elastic cc	€ 168,94

SatinSheet DrawSheet

Combine SatinSheet DrawSheet with SatinSheet 2Direction BaseSheet for low-friction side to side

IM4116S	SatinSheet DrawSheet - mini	W2000xL1000 mm. Stripes	€ 56,77
IM4118S	SatinSheet DrawSheet - midi	W2000xL1400 mm. Stripes	€ 59,07

SatinSheet 4Direction

Combine SatinSheet 4Direction with SatinSheet 2Direction Basesheet for low-friction in all four directions

IM4115S	SatinSheet 4Direction - midi	W2000xL1400 mm (satin W1470xL1400 mm). Check	€ 114,66
IM4115Sf	SatinSheet 4Direction w/handles - midi	W2000xL1400 mm (satin W1470xL1400 mm). Check	€ 149,70
IM4121S	SatinSheet 4Direction - maxi	W2000xL2000 mm (satin W1470xL2000 mm). Check	€ 180,00
IM4121Sf	SatinSheet 4Direction w/handles - maxi	W2000xL2000 mm (satin W1470xL2000 mm). Check	€ 168,94

SatinSheet In2Sheet 4Direction

Combine SatinSheet In2Sheet 4Direction with SatinSheet 2Direction Basesheet for low-friction in all four directions

IM4110S	In2Sheet 4Direction w/wings - mini	W850xL900 mm. Back side: Satin. Top side: Micro. Wing	€ 102,03
IM4110Sf	In2Sheet 4Direction w/wings + handles - n	W850xL900 mm. Back side: Satin. Top side: Micro. Wing	€ 111,76
IM4140S	In2Sheet 4Direction w/wings - midi	W1200xL1400mm. Back side: Satin. Top side: Micro. Wir	€ 161,61
IM4140Sf	In2Sheet 4Direction w/wings + handles - n	W1200xL1400mm. Back side: Satin. Top side: Micro. Wir	€ 206,09

SatinSheet Pillow Case

IM4117S	SatinSheet 2Direction PillowCase	W700xL600 mm (satin W700xL600 mm). Stripes	€ 44,53
---------	----------------------------------	--	---------

In2Sheet

IM4109	In2Sheet w/handles - mini	W850xL900 mm. Back side: Satin. Top side: Micro	€ 76,76
IM4128	In2Sheet w/handles - midi	W1050xL1250 mm. Back side: Satin. Top side: Micro	€ 103,29

4WayGlide LPL System

4WayGlide Mattress

IM140/20C	4WayGlide Mattress LPL	W1400xL2100 mm. Top side: Polycotton	€ 234,35
IM140/20C	4WayGlide Mattress LPL	W1400xL2100 mm. Top side: PU	€ 309,81

4WayGlide NylonSheet

IM85/200 4WayGlide NylonSheet	W850xL2000 mm	€ 131,55
IM85/200 4WayGlide NylonSheet	W850xL2000 mm. Non-Slip sides	€ 142,13
IM90/200 4WayGlide NylonSheet	W900xL2050 mm. Non-Slip sides	€ 114,51
IM85/200L 4WayGlide NylonSheet LPL w/lock	W850xL2000 mm.	€ 170,46
IM85/200L 4WayGlide NylonSheet LPL w/lock	W850xL2000 mm. Non-Slip sides	€ 215,96
IM85/200 4WayGlide NylonSheet LPL w/2x lock	W850xL2000 mm. Non-Slip sides	€ 258,36
IM85/200 4WayGlide NylonSheet LPL w/3x lock	W850xL2000 mm. Non-Slip sides	€ 302,50
<i>Complete system Includes 1pc of each Mattress and NylonSheet</i>		

4WayGlide LPL Cover

IM140/200 4WayGlide Cover	W1400xL2200 mm. Top side: cotton. Back side: PU	€ 176,97
---------------------------	---	----------

4WayGlide Mattress Maxi

IM180/200 4WayGlide Mattress LPL - maxi	W1800xL2000 mm. Top side: Polycotton	€ 382,65
IM180/200 4WayGlide Mattress LPL - maxi	W1800xL2000 mm. Top side: PU	€ 519,77

4WayGlide NylonSheet Maxi

IM110/200 4WayGlide NylonSheet LPL w/lock - maxi	W1100xL2000 mm. Non-Slip sides	€ 223,32
IM120/200 4WayGlide NylonSheet	W1200xL2050 mm.	€ 79,38
IM120/200 4WayGlide NylonSheet LPL w/3x lock - m	W1200xL2000 mm. Non-Slip sides	€ 327,72

Complete system includes 1pc of Mattress and NylonSheet

TwinSheet4Glide

IM150/200 TwinSheet4Glide - long	W1500xL2000 mm. Absorbing pad	€ 373,99
IM150/140 TwinSheet4Glide - down	W1500xL1400 mm. Absorbing pad	€ 249,07
IM150/75 TwinSheet4Glide - top	W1500xL850 mm. Absorbing pad	€ 189,68
IMLWS60 Loop Washing Strap	W60xL520 mm	€ 14,52

Complete system includes 1pc of Mattress and NylonSheet

4WayGlide Mini

IM140/100 4WayGlide Mattress - mini	W1400x1000 mm. Top side: Polycotton	€ 106,24
IM140/100 4WayGlide NylonSheet - mini	W1600xL1000 mm. Non-Slip sides	€ 101,93
<i>Complete system includes 1pc of Mattress and NylonSheet</i>		

4WayGlide Midi

IM140/150 4WayGlide Mattress w/lock - midi	W1400x1500 mm. Top side: Polycotton	€ 170,22
IM140/150 4WayGlide NylonSheet - midi	W1600xL1500 mm. Non-Slip sides	€ 180,73
<i>Complete system includes 1pc of Mattress and NylonSheet</i>		

MultiGlide high-low friction (open on the width)

(width x length)

IM100	MultiGlide vertical XXS	W1000xL700 mm.	€ 55,66
IM103/1	MultiGlide vertical XS	W1000xL1350 mm	€ 65,82
IM102	MultiGlide vertical S	W1000xL200 mm	€ 103,43
IM108	MultiGlide vertical M	W1080xL800 mm	€ 53,87
IM120	MultiGlide vertical L	W1200xL700 mm	€ 60,50
IM130	MultiGlide vertical XXL	W1300xL900 mm	€ 63,89
IM168	MultiGlide vertical XXXL	W1600xL800 mm	€ 102,46
IM200	MultiGlide vertical XXXXL	W2000xL700 mm	€ 80,95
IM40	MultiGlide lateral XXS	W400xL400 mm	€ 32,28
IM62	MultiGlide lateral XS	W600xL250 mm	€ 39,54
IM66	MultiGlide lateral S	W600xL600 mm.	€ 35,21
IM77	MultiGlide lateral M	W720xL700 mm	€ 37,12
IM88	MultiGlide lateral L	W800xL800 mm	€ 87,94
IM80	MultiGlide lateral XL	W720xL1000 mm	€ 41,33
IM83	MultiGlide lateral XXL	W720xL1300 mm	€ 66,86

IM82	MultiGlide lateral XXXL	W700xL2000 mm	€ 67,49
IM100/140	MultiGlide open	W1000xL1400 mm	€ 41,62

MultiGlide SG (SpilerGlide) low-low friction (open on the w (width x length)

IM00100	MultiGlide SG vertical XXS	W1000xL700 mm	€ 51,64
IM00108	MultiGlide SG vertical M	W1080xL800 mm	€ 71,00
IM00120/1	MultiGlide SG vertical XL	W1200xL900 mm	€ 71,34
IM00139	MultiGlide SG vertical XXL	W1300xL900 mm	€ 92,78
IM00130	MultiGlide SG vertical XXL - mini	W1300xL800 mm	€ 73,86
IM00170	MultiGlide SG vertical XXL - wide	W1400xL1200 mm	€ 95,20
IM00200	MultiGlide SG vertical XXXXL	W2000xL700 mm	€ 105,37
IM0040	MultiGlide SG lateral XXS	W400xL400 mm	€ 30,13
IM0066	MultiGlide SG lateral S	W600xL600 mm	€ 34,70
IM0077	MultiGlide SG lateral M	W720xL700 mm	€ 37,12
IM0088	MultiGlide SG lateral L	W800xL800 mm	€ 48,06
IM0080	MultiGlide SG lateral XL	W720xL1000 mm	€ 46,95
IM0083	MultiGlide SG lateral XXL	W720xL1300 mm	€ 47,34

MultiGlide with handles

IM130DR	MultiGlide with handles	W1300xL800 mm	€ 91,55
IM160DR	MultiGlide with handles	W1600xL800 mm	€ 101,88

Single Patient Use

MultiGlide Single Patient Use

IM301	MultiGlide Single Patient Use lateral XXL	W730xL1450 mm. 5 boxes of 10 pcs	€ 1.072,01
IM302	MultiGlide Single Patient Use vertical XXL	W1300xL900 mm. 5 boxes of 10 pcs	€ 1.721,15
IM303	MultiGlide Single Patient Use vertical L	W1200xL730 mm. 5 boxes of 10 pcs	€ 911,95
IM304	MultiGlide Single Patient Use open L	W2000xL850 mm. 5 boxes of 20 pcs	€ 2.614,47
IM305	MultiGlide Single Patient Use open M	W1450xL1000 mm. 5 boxes of 20 pcs	€ 1.500,79
IM306	MultiGlide Single Patient Use open S (2 pc	W1000xL1000 mm. 5 boxes of 10 pcs	€ 1.460,13

Glide Cushion & TurningMattress

IM60/50	Glide Cushion 600x500 mm	nylon inside top side cotton green	€ 85,52
IM60/90	GlideCushion 600x900 mm	nylon inside top side cotton green	€ 101,83
IM69/70	GlideCushion 690x700 mm	nylon inside top side cotton green	€ 92,30
IM90/90	GlideCushion 900x900 mm	nylon inside top side cotton green	€ 147,77
IM92/195	TurningMattress 920x1950 mm	nylon inside top side cotton green	€ 337,20

Glide Cushion Set

IM4303	GlideCushion X-Small Set	incl IM60/50 + IM6000N	€ 126,57
--------	--------------------------	------------------------	----------

Positioning - Sitting

Swivel cushions

IM99999	Autoturn	soft auto cushion w belt NS-bottom velour top, 450 mm d	€ 138,98
IM999405	EasyTurn	soft swivel cushion NS-bottom velour top , 400 mm dimat	€ 91,55
IM990454	EasyTurn	soft swivel cushion NS-bottom velour top , 450 mm dimat	€ 92,52
IM999505	EasyTurn	soft swivel cushion NS-bottom velour top , 500 mm dimat	€ 100,14

OneWayGlide

IM39/8	OneWayGlide NS	370x430 mm tube w non-slip top grey	€ 71,75
IM38/8	OneWayGlide NS	450x450 mm tube w non-slip top grey	€ 73,11
IM45/8	OneWayGlide Velour	370x430 mm tube velour top	€ 93,53
IM50/8	OneWayGlide Velour	500x400 mm tube velour top	€ 100,26
IM44/8	OneWayGlide Velour Pad	370x430 mm one layer pad velour top	€ 45,93
IM33/8	OneWayGlide Long w Straps	400x1300 mm long open model velour top	€ 99,97
IM34/8	OneWayGlide Long	400x1300 mm long open model velour top	€ 105,08

IM51/8	OneWayGlide Long	450x1300 mm long open model velour top	€ 107,79
IM35/8	OneWayGlide Long	550x1300 mm long open model velour top	€ 130,08

Transfer Mattress

IM60/195	Transfer Mattress Narrow	nylon inside top side cotton green, 600x1950 mm	€ 176,71
----------	--------------------------	---	----------

Transfer Mattress cover

IM60220N	Transfer Mattress Cover	600x2200 mm waterproof nylon (for IM60/195) green	€ 67,42
----------	-------------------------	---	---------

Transfer Mattress Set

IM4304	Transfer Mattress Narrow Set	incl IM60/195 + IM60220N	€ 250,28
--------	------------------------------	--------------------------	----------

Support Sling

IM428	Sling - long	W160xL1900 mm. W/2x 4 handles	€ 80,15
-------	--------------	-------------------------------	---------

SupportBelt

IM800900	SupportBelt XXS	Size: 450-900 mm. 4 handles. Pink label	€ 61,95
IM801050	SupportBelt XS	Size: 650-1050 mm. 4 handles. Light Blue label	€ 65,29
IM801051	SupportBelt S	Size: 650-1050 mm. 4 handles. Red label	€ 71,15
IM801301	SupportBelt M	Size: 900-1300 mm. 6 handles. Yellow label	€ 80,05
IM801501	SupportBelt L	Size: 1200-1500 mm. 8 handles. Green label	€ 93,22
IM801701	SupportBelt XL	Size: 1300-1700 mm. 11 handles. Blue label	€ 110,79

SmartMove

IM426	SmartMove w/4 handles	W210xL600 mm	€ 53,87
-------	-----------------------	--------------	---------

BedString

IM429	BedString w/handles	W85xL1700 mm. Loop fixing	€ 66,70
-------	---------------------	---------------------------	---------

Special Solutions

SlingOn

MiniBoard

IM801505 MiniBoard L	W540xL1500 mm. W/6 handles	€ 225,88
IM801245 MiniBoard w/head support L	W540xL1240 mm. W/7 handles	€ 208,85

OneManSling

IM460PU OneManSling	W75xL3000 mm. Sling w/PU cushion	€ 105,10
---------------------	----------------------------------	----------

Positioning Wedge

IM53NS Positioning Wedge S	W300xL530xH200 mm. Incl. Non-Slip cover	€ 111,27
IM53/80N Positioning Wedge L		€ 157,30
IM5300NS Positioning Wedge Cover	W300xL530xH200 mm. For IM53NS	€ 67,23

MultiGlide Glove

IM35 MultiGlide Glove	W200xL500 mm. 1 pair. Nylon	€ 15,88
-----------------------	-----------------------------	---------

Non-Slip

IM2248 Non-Slip Square	W380xL480xH5 mm	€ 27,93
IM2246 Non-Slip Roll. 5 metres	W1500xL5000 mm	€ 165,72